



**SEARCH MEDICAL Co.,Ltd.**


16, Cheomdan venture so-ro 37beon-gil,  
Buk-gu, Gwangju, Republic of Korea [61003]  
Tel. 062-714-2699 Fax. 070-4855-2609

# PDO

POLYDIOXANONE



SEARCH MEDICAL  
(주)서치메디컬



# PDO

## POLYDIOXANONE

It is a PDO suture made of Polydioxanone and is an absorbable suture. To suture the tissue, it is attached to a stainless steel needle that can pass through the tissue and is absorbed after a period of use.


PDO suture is a single structure with a smooth surface and good flexibility, so it is easy to handle and has less resistance when passing through the tissue.



# Polydioxanone

Absorbable Suture

Polydioxanone Violet, Mono filament



1.Raw material  
Polydioxanone

5.Validity  
3years from date of manufacture

2.Color  
Violet

3.Structure  
Mono filament

4.Suture with Needle



[USP]

7/0, 6/0, 5/0, 4/0, 3/0, 2/0, 0,1, 2

[Length]

70cm

Needle Identification

TR= Taper Point (Round)

BODY	POINT

**TR**

CC=Conventional cutting

BODY	POINT

**CC**

RC=Reverse cutting

BODY	POINT

**RC**

TB= Taper Point (Blunt)

BODY	POINT

**TB**

SC=Spatula Cutting

BODY	POINT

**SC**




MS=Micro point spatula














BODY	POINT


















**MS**

SEARCH MEDICAL























Needle Shape [Circle]

























Needle Type	Suture	Needle Product		
	USP	Needle Length [mm]	Circle	Model Name
 	 7-0 Violet 700mm	Length [mm]		
 		8	1/2	STR7008
 		8	1/2	STB7008
 		8	1/2	SRC7008
 		8	3/8	SSC7008
 		8	3/8	SMS7008
































Needle Type	Suture	Needle Product		
	USP	Needle Length [mm]	Circle	Model Name
 	 6-0 Violet 700mm	Length [mm]		
 		8	1/2	STR6008
 		8	1/2	STB6008
 		8	1/2	SRC6008
 		8	3/8	SSC6008
 		8	3/8	SMS6008
 		10	3/8	SRC6010
 		12	3/8	SRC6012

SEARCH MEDICAL














Needle Type	Suture	Needle Product		
	USP	Needle Length [mm]	Circle	Model Name
 	5-0 Violet 700mm	Length [mm]		
 		15	1/2	STR5015
 		15	1/2	STB5015
 		18	1/2	STR5018
 		18	1/2	STB5018
 		10	3/8	SRC5010
 		12	3/8	SRC5012
 		13	3/8	STB5013
 		16	3/8	SRC5016
 		16	3/8	STR5016
 	17	3/8	STB5017	




























Needle Type	Suture	Needle Product		
	USP	Needle Length [mm]	Circle	Model Name
 	4-0 Violet 700mm	Length [mm]		
 		18	1/2	STB4018
 		20	1/2	STR4020
 		20	1/2	STB4020
 		22	1/2	STR4022
 		22	1/2	STB4022
 		12	3/8	SRC4012
 		16	3/8	SRC4016
 		16	3/8	STR4016
 		16	3/8	STB4016
 	19	3/8	SRC4019	
























Needle Type	Suture	Needle Product			
	USP	Needle Length [mm]	Circle	Model Name	
 	 3-0 Violet 700mm	Length [mm]	15	1/2	STR3015
 			15	1/2	STB3015
 			18	1/2	STR3018
 			18	1/2	STB3018
 			20	1/2	STR3020
 			20	1/2	STB3020
 			22	1/2	STR3022
 			22	1/2	STB3022
 			25	1/2	SSC3025
 			26	1/2	STR3026
 			26	1/2	STB3026
 			30	1/2	SCC3030
 			30	1/2	STR3030
 			35	1/2	STR3035
 			35	1/2	STB3035















SEARCH MEDICAL
























Needle Type	Suture	Needle Product			
	USP	Needle Length [mm]	Circle	Model Name	
 	 3-0 Violet 700mm	Length [mm]	16	3/8	SRC3016
 		16	3/8	STR3016	
 		16	3/8	STB3016	
 		24	3/8	SRC3024	
 		30	3/8	SRC3030	
 		45	Straight	SCC3045	

Needle Type	Suture	Needle Product			
	USP	Needle Length [mm]	Circle	Model Name	
 	 2-0 Violet 700mm	Length [mm]	18	1/2	STR2018
 		18	1/2	STB2018	
 		20	1/2	STR2020	
 		20	1/2	STB2020	
 		22	1/2	STR2022	
 		22	1/2	STB2022	
 		26	1/2	STR2026	
 		26	1/2	STB2026	
 		30	1/2	SCC2030	
 		35	1/2	STR2035	
 		35	1/2	STB2035	
 		50	1/2	STR2050	
 		50	1/2	STB2050	























Needle Type	Suture	Needle Product		
	USP	Needle Length [mm]	Circle	Model Name
 	 2-0 Violet 700mm	Length [mm]		
 		16	3/8	STR2016
 		16	3/8	STB2016
 		16	3/8	SCC2016
 		19	3/8	SRC2019
 		24	3/8	SRC2024
 		28	3/8	SCC2028
 		30	3/8	SRC2030
 		50	3/8	SCC2050
 		45	Straight	SCC2045



















SEARCH MEDICAL

Needle Type	Suture	Needle Product		
	USP Length [mm]	Needle Length [mm]	Circle	Model Name
	 Violet 700mm	20	1/2	STR020
		20	1/2	STB020
		22	1/2	STR022
		22	1/2	STB022
		26	1/2	STR026
		26	1/2	STB026
		30	1/2	SCC030
		35	1/2	STR035
		35	1/2	STB035
		40	1/2	STR040
		40	1/2	STB040
		50	1/2	STR050
		50	1/2	STB050

Needle Type	Suture		Needle Product	
	USP	Length [mm]	Circle	Model Name
 	 Violet 700mm	24	3/8	SRC024
 		28	3/8	SCC028
 		30	3/8	SRC030
 		39	3/8	STR039
 		39	3/8	STB039
 		40	3/8	SCC040
 		45	3/8	STR045
 		45	3/8	STB045
 		50	3/8	SCC050
 		45	Straight	SCC045
 		60	Straight	SCC060

SEARCH MEDICAL

Needle Type	Suture		Needle Product		
	USP	Length [mm]	Needle Length [mm]	Circle	Model Name
 	1 Violet 700mm	Length [mm]	26	1/2	STR126
 			26	1/2	STB126
 			30	1/2	SCC130
 			35	1/2	STR135
 			35	1/2	STB135
 			40	1/2	STR140
 			24	3/8	SRC124
 			28	3/8	SCC128
 			30	3/8	SRC130
 			45	Straight	SCC145
 	60	Straight	SCC160		

Needle Type	Suture	Needle Product				
	USP	Needle Length [mm]	Circle	Model Name		
 	 Violet 700mm	Length [mm]	40	1/2	STB240	
			50	1/2	STR250	
			50	1/2	STB250	
 		 Violet 700mm	Length [mm]	39	3/8	STR239
				 	39	3/8
 		 Violet 700mm	Length [mm]	40	3/8	SCC240
 				45	3/8	STR245
		 	 Violet 700mm	Length [mm]	45	3/8
 					50	3/8

SEARCH MEDICAL